



ARJO SARA 3000 TRANSFER SLING

\$801.00



SKU: TSS.511-L1

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Arjo universal transfer sling is designed for use with specific standing and raising aids. Universal transfer sling is intended not only to provide support for standing safely but also to promote active participation by the patient in standing movement. Universal transfer sling is easy to use while transferring the resident from the bed to a wheelchair. It comes in one size with a Safe Working Load of 200 kilograms.

FEATURES

- Designed for use with specific standing and raising aids.
- Universal transfer sling is intended not only to provide support for standing safely but also to promote active participation by the patient in standing movement.
- Easy to use while transferring the resident from the bed to a wheelchair.
- One size fits all

SPECIFICATIONS

- Safe working load: 200 kg/440 lbs