



ARJO SARA FLEX - SIT TO STAND AID



SKU: HEB0000-22-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Arjo Sara Flex is specifically designed to simplify everyday transfer and care tasks. With the Sara Flex, a single caregiver can effortlessly assist a patient or resident from a seated to a standing position in one smooth and natural movement.

The Arjo Sara Flex is available in two versions: with or without a scale. It comes equipped with two batteries and a charger, ensuring continuous and reliable operation. With a maximum weight capacity of 200kg, it provides robust support for various individuals.

Patient Support: The Sara Flex prioritizes the comfort and security of the patient throughout the sit-to-stand movement. Its design promotes a sense of ease and stability, ensuring a reassuring experience for the patient.

Facilitating a Range of Key Transfers and Care Tasks: With the Sara Flex, performing essential transfers and care tasks becomes effortless for the single caregiver. This innovative device requires no adjustments, enabling efficient and seamless operations.

Natural Forward Movement for the Patient: The Arjo Sara Flex facilitates a natural and comfortable sit-to-stand movement through its forward motion feature. Combined with the flexible leg support, this design ensures a smooth and ergonomic experience for the patient.

DSDC DEMENTIA DESIGN ACCREDITED PRODUCT

The Dementia Services Development Centre (DSDC) at the University of Stirling has reviewed Sara Flex and awarded it a DSDC Product Accreditation Rating Class 1B.

FEATURES

- **Flexible Silicone Leg Support:** The Sara Flex features a curved design with a soft and flexible surface that supports the lower legs and knees throughout the sit-to-stand movement. This large surface requires no adjustments and caters to a broad range of patients while delivering good comfort.
- **Flexible Silicone Leg Strap:** Positioned low and offering multiple adjustment options, the Arjo Sara Flex ensures that patients can be securely positioned from a seated to a standing position.
- **Easy-to-Use Hand Control with Display:** The hand control of the Arjo Sara Flex is conveniently located on a long flexible cord, allowing the caregiver to stay close to the patient throughout the transfer. It also features a display for easy monitoring and control.
- **Caregiver Handle:** The Arjo Sara Flex is equipped with a curved caregiver handle that enables the caregiver to choose the most comfortable and convenient holding position.
- **Patient Handle:** Featuring an angled and contrast-coloured design, the patient handle of the Arjo Sara Flex positions the hands, wrists, and elbows of the patient optimally during the sit-to-stand movement.
- **Low Horizontal Footplate:** The Arjo Sara Flex is designed with a low horizontal footplate, keeping the patient's feet close to the ground. This design prevents the knees from being raised too high, making it easier for the patient to stand.
- **Scale (Option):** The Arjo Sara Flex offers an optional built-in electronic scale that shows the patient's weight on the display of the hand control, providing convenient weight measurement during transfers.

SPECIFICATIONS

Length	960mm
Height	1160mm
External Width Legs Open	1030mm
External Width Legs Closed	660mm
Internal Width Legs Open	900mm
Internal Width Legs Closed	530mm
SWL	200kg

VARIATIONS AND ACCESSORIES

Sara Flex without Scale - comes with 2 batteries and 1 charger (SWL 200kg).

Sara Flex with Scale - comes with 2 batteries and 1 charger (SWL 200kg).

BROCHURES

[Sara flex Brochure](#)

[Sara Flex Specification Sheet](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au