




ARJO SARA PLUS - SIT TO STAND AID

\$18,781.00 - \$23,621.00 Price range: **\$18,781.00 through \$23,621.00**



SKU: HEP0001-AU-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS OF THE ARJO SARA PLUS

The Arjo Sara Plus is an ergonomic standing and raising aid that mobilises patients during everyday activities. The Sara Plus can be operated by a single caregiver for transfers and patient training. Equipped with a comprehensive range of safety features, this device drastically reduces the risk of injuries associated with manual patient handling, ensuring a safer and more dignified experience for both the patient and caregiver.

Additionally, the Arjo Sara Plus offers advanced features and customizable options to accommodate the unique needs of each patient. Its intuitive control panel allows caregivers to adjust the lifting speed, height, and support levels, ensuring optimal comfort and stability during transfers. The device is designed with patient dignity in mind, providing a smooth and gentle lifting motion that minimizes discomfort and maintains a natural posture. With its versatility and user-friendly design, the Arjo Sara Plus promotes independence and enhances the overall quality of care for individuals with limited mobility.

FEATURES

- The Arc-Rest provides upper body support during the raising action
- Designed to offer balance, stepping and walking training
- The first-step system, incorporating the Proactive Pad™ and a detachable footplate, enables a range of training exercises directly on the floor
- Dual controls on the mast and the hand control. The hand control enables the caregiver to move and stay close to the patient and maintain good eye contact throughout the process.
- A detachable footplate enables stepping and gait training with the patient's feet directly on the floor.

- Comes with 2 batteries and 1 charger
- Optional integrated scale contributes to an easier weighing routine

VIDEO

[embed]https://www.youtube.com/watch?v=00iKQKDH_2o[/embed]

SPECIFICATIONS

- Max. safe working load - 420 lbs (190 kg)
- Weight - 187 lbs (85 kg)
- Weight with Scale - 209 lbs (95 kg)
- Total length - 40 5/8" (1033 mm)
- Min. lifting height - 35 5/8" (906 mm)
- Max. lifting height - 58 7/8" (1497 mm)
- Total height - 47" (1195 mm)
- Height of chassis - 4 5/8" (116 mm)
- External width legs closed - 25 3/8" (646 mm)
- Internal width legs open - 39 1/2" (1002 mm)
- Low friction castors, the rear two with brakes
- Battery 24V 5Ah
- Two batteries and one charger are supplied with each Sara Plus
- Protection class - Lifter IPX4
- Handset IP67

BROCHURES

[Sara Plus Specification Sheet](#)

[Arjo Sara Plus Cord Disinfection Clinical Summary](#)