

**ARJO SARA STEDY - SIT TO STAND AID****\$3,388.00****SKU: NTB2000-H****PRODUCT DESCRIPTION****PURPOSE AND BENEFITS**

Arjo Sara Stedy and Arjo Sara Stedy Compact enable a single caregiver to assist patients or residents in performing stand transfers throughout the day. With their ergonomic design and sturdy construction, Sara Stedy and Sara Stedy Compact are essential tools for healthcare facilities aiming to enhance patient mobility and reduce the risk of injuries. These innovative transfer aids provide efficient and safe sit-to-stand transitions, promoting better patient care and ensuring a positive experience for both caregivers and patients.

**FEATURES**

- Arjo Sara Stedy and Sara Stedy Compact are highly effective in promoting early mobilization by allowing patients to engage in sit-to-stand transfers at an early stage of their recovery or rehabilitation process.
- With the option of two sizes, Standard and Compact, healthcare facilities can choose the most suitable model based on their specific spatial requirements and patient needs, ensuring optimal comfort and usability.
- Being a non-powered device, Sara Stedy and Sara Stedy Compact provide the flexibility of use at any time, without the need for electrical outlets or charging. This feature makes them reliable and readily available for caregivers to assist patients with sit-to-stand transfers whenever necessary, promoting efficiency and convenience in healthcare settings.

**DSDC DEMENTIA DESIGN ACCREDITED PRODUCT**

The [Dementia Services Development Centre](#) at the University of Stirling has reviewed Sara Stedy and awarded it a DSDC Product Accreditation Rating Class 1B.

## VIDEO

[embed]<https://www.youtube.com/watch?v=iRlsA7w7o2g>[/embed]

## SPECIFICATIONS

### SARA STEDY

Designed for patients 149-193 cm (4'11"-6'6") tall and weighing up to 182 kg (400 lb).

<b>Length</b>	920mm
<b>Height</b>	1050mm
<b>External Width   Legs Open</b>	880mm
<b>External Width   Legs Closed</b>	490mm
<b>Internal Width   Legs Open</b>	750mm
<b>Internal Width   Legs Closed</b>	360mm
<b>User Height</b>	149cm - 193cm
<b>SWL</b>	182kg

### SARA STEDY COMPACT

Designed for smaller patients, 138-172 cm (4'6"-5'8") tall and weighing up to 182 kg (400 lb).

<b>Length</b>	890mm
<b>Height</b>	960mm
<b>External Width   Legs Open</b>	880mm
<b>External Width   Legs Closed</b>	490mm
<b>Internal Width   Legs Open</b>	750mm
<b>Internal Width   Legs Closed</b>	360mm
<b>User Height</b>	138cm - 172cm
<b>SWL</b>	182kg

Buy this Arjo Sara Steady online at Medimart using Stripe, PayPal, Afterpay or [NDIS method](#) today!