




ARJO SARA STEDY - SIT TO STAND AID



SKU: NTB2000-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Arjo Sara Stedy and Arjo Sara Stedy Compact enable a single caregiver to assist patients or residents in performing stand transfers throughout the day. With their ergonomic design and sturdy construction, Sara Stedy and Sara Stedy Compact are essential tools for healthcare facilities aiming to enhance patient mobility and reduce the risk of injuries. These innovative transfer aids provide efficient and safe sit-to-stand transitions, promoting better patient care and ensuring a positive experience for both caregivers and patients.

FEATURES

- Arjo Sara Stedy and Sara Stedy Compact are highly effective in promoting early mobilization by allowing patients to engage in sit-to-stand transfers at an early stage of their recovery or rehabilitation process.
- With the option of two sizes, Standard and Compact, healthcare facilities can choose the most suitable model based on their specific spatial requirements and patient needs, ensuring optimal comfort and usability.
- Being a non-powered device, Sara Stedy and Sara Stedy Compact provide the flexibility of use at any time, without the need for electrical outlets or charging. This feature makes them reliable and readily available for caregivers to assist patients with sit-to-stand transfers whenever necessary, promoting efficiency and convenience in healthcare settings.

DSDC DEMENTIA DESIGN ACCREDITED PRODUCT

The [Dementia Services Development Centre](#) at the University of Stirling has reviewed Sara Stedy and awarded it a DSDC Product Accreditation Rating Class 1B.

VIDEO

[embed]<https://www.youtube.com/watch?v=iRlsA7w7o2g>[/embed]

SPECIFICATIONS

SARA STEDY

Designed for patients 149-193 cm (4'11"-6'6") tall and weighing up to 182 kg (400 lb).

Length	920mm
Height	1050mm
External Width Legs Open	880mm
External Width Legs Closed	490mm
Internal Width Legs Open	750mm
Internal Width Legs Closed	360mm
User Height	149cm - 193cm
SWL	182kg

SARA STEDY COMPACT

Designed for smaller patients, 138-172 cm (4'6"-5'8") tall and weighing up to 182 kg (400 lb).

Length	890mm
Height	960mm
External Width Legs Open	880mm
External Width Legs Closed	490mm
Internal Width Legs Open	750mm
Internal Width Legs Closed	360mm
User Height	138cm - 172cm
SWL	182kg

Buy this Arjo Sara Steady online at Medimart using Stripe, PayPal, Afterpay or [NDIS method](#) today!

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au