



ARJO SORRENTO - FLOOR BED



SKU: LAU-G420-88

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Sorrento Floor Bed is a stylish, yet economical choice for long-term care, aged care, and all nursing and community healthcare environments. Designed to reduce the risk of harm from falls whilst allowing you to provide a premium aged care bed for your residents and peace of mind for their families and caregivers. With high-end customisation options and features available, the Sorrento Floor Bed suits any environment.

FEATURES

- Electrically operated backrest and height adjustment
- 4-section bed | backrest and knee-break
- Mechanical calf adjustment | 4-section model only
- Battery back-up
- Fitted base extension
- Fully moulded undercarriage cover
- Removable sheet metal platform sections
- Easy to clean ABS moulded head and foot panels
- 125mm central locking system – 2 x brake | 1 x track lock | 1 plain
- Available in Single and King Single Sizes

SPECIFICATIONS

Sorrento Floor Bed Single

Overall Length	2115 mm
Overall Length (base extension extended)	2310 mm
Overall Width	9200 mm
Overall Width (rails fitted)	1025 mm

Mattress Surface Length (between panels)	2070 mm
Bed Platform Height Range	1700 mm – 690 mm (+/- 5 mm)
Back Rest Angle Adjustment	0° to 70°
Knee Bend Angle Adjustment	0° to 35°
Trendelenberg Tilt Adjustment (head and foot down)	0° to 14°
Safe Working Load	250 kg

Sorrteno Floor Bed King Single

Overall Length	2115 mm
Overall Length (base extension extended)	2310 mm
Overall Width	1090 mm
Overall Width (rails fitted)	1200 mm
Mattress Surface Length (between panels)	2070 mm
Bed Platform Height Range	170 mm to 690 mm (+/- 5mm)
Back Rest Angle Adjustment	0° to 70°
Knee Bend Angle Adjustment	0° to 35°
Trendelenberg Tilt Adjustment (head and foot down)	0° to 14°
Safe Working Load	250 kg