



ASPIRE 600MM WEIGHT BEARING PLATFORM SHOWER COMMODE

\$2,044.00 Original price was: \$2,044.00.<u>\$1,840.00</u>Current price is: \$1,840.00.



SKU: BTC061400

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aspire Weight Bearing Platform Shower Commode is a high-quality mobility aid designed to make bathing and toileting easier for those with limited mobility. With a weight capacity of up to 400 kg, this commode is perfect for bariatric patients who require additional support and stability.

The platform design provides extra space and stability for users to comfortably sit while showering or using the commode. The durable, stainless steel frame is rust-resistant and designed to withstand daily use. The commode features a padded seat and backrest for added comfort and support, as well as swing-away footrests for easy transfer in and out of the chair.

The Aspire 600mm Weight Bearing Platform Shower Commode is easy to clean and maintain, making it a practical and hygienic solution for any care setting.

FEATURES

- Heavy duty reinforced wide frame can support bariatric users up to 400 kg
- 4 x total locking stainless steel castors provide brake access from all angles
- Height adjustable reinforced armrests offer drop flat
- <u>Seat</u>, pan, and pan carrier not included sold separately

SPECIFICATIONS

Code Display Name

Seat Width Safe Working Load



Code Display Name Seat Width Safe Working Load

BTC061400 Attendant Propelled Commode - Weight Bearing 600 mm 400 kg