




ASPIRE 600MM WEIGHT BEARING PLATFORM SHOWER COMMODE



SKU: BTC061400

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aspire Weight Bearing Platform Shower Commode is a high-quality mobility aid designed to make bathing and toileting easier for those with limited mobility. With a weight capacity of up to 400 kg, this commode is perfect for bariatric patients who require additional support and stability.

The platform design provides extra space and stability for users to comfortably sit while showering or using the commode. The durable, stainless steel frame is rust-resistant and designed to withstand daily use. The commode features a padded seat and backrest for added comfort and support, as well as swing-away footrests for easy transfer in and out of the chair.

The Aspire 600mm Weight Bearing Platform Shower Commode is easy to clean and maintain, making it a practical and hygienic solution for any care setting.

FEATURES

- Heavy duty reinforced wide frame can support bariatric users up to 400 kg
- 4 x total locking stainless steel castors provide brake access from all angles
- Height adjustable reinforced armrests offer drop flat
- **[Seat](#), pan, and pan carrier not included - sold separately**

SPECIFICATIONS

Code

Display Name

Seat Width Safe Working Load

Code	Display Name	Seat Width	Safe Working Load
BTC061400	Attendant Propelled Commode - Weight Bearing Platform - 600 mm Wide	600 mm	400 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au