



ASPIRE A205 PATIENT LIFTER



SKU: LSS390610

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aidacare Aspire A205 Patient Lifter boasts a compact yet versatile design, making it a top choice for diverse settings ranging from homes to community facilities. Crafted from durable aluminium, this mobile lifter comes equipped with cutting-edge technology and features that prioritise both patient and caregiver safety.

FEATURES

- **Easy Storage and Transportation:** Quickly disassemble the lifter without needing any tools, making storage and transportation a breeze.
- **Ergonomic Design:** The lifter's Linak slimline technology-based handset is user-friendly, featuring an intuitive control panel. Moreover, built-in battery and service indicator lights aid in timely maintenance.
- **Flexibility in Usage:** Caregivers can choose from three height settings, ensuring that the lifter caters to different patient needs and scenarios.
- **Enhanced Stability:** The powered leg function enables caregivers to open and close the lifter's legs, ensuring stability during lifts.
- **Swift Yoke Attachment:** Designed with busy caregivers in mind, the one-touch quick-release yoke attachment also features padding for the patient's safety.
- **Durable Construction:** This lifter stands out with its strong powder-coated aluminium frame and industry-leading Linak electronics, offering a robust warranty provision.
- **Safety First:** Pivotal safety features include an emergency stop button, antientrapment technology, and emergency lowering. Safety is paramount for both the caregiver and the patient.

- **Advanced Diagnostics:** The state-of-the-art Linak Jumbo Care control system tracks the lifter's usage and service requirements, ensuring compliance and safety.
- **Optimal Mobility:** The Tente casters comprise 2 front castors (75 mm) for optimal under-bed clearance and 2 rear locking castors (100 mm) to enhance safety and manoeuvrability.
- **Compliance:** This lifter complies with AS/NZS ISO 10535:2011, the established Australian standard for hoists.

The Aspire A205 Patient Lifter is a harmonious blend of efficiency, safety, and technological prowess. Designed to cater to the needs of both caregivers and patients, this lifter ensures a smooth and secure lifting experience. Its compact nature combined with its myriad of features makes it an indispensable asset in varied care settings.

SPECIFICATIONS

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Lift Range	Leg Gait Range	Electric Leg
LSS390610	Aspire A205 Patient Lifter	1135 mm	620 - 1095 mm	2045 mm	205 kg	Aluminium	38.1 kg	395 - 1805 mm	547 - 1015 mm	True

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au