

ASPIRE ACTIVE AIR 8 - ALTERNATING MATTRESS REPLACEMENT

\$3,718.00 - \$4,864.00 Price range: \$3,718.00 through \$4,864.00



SKU: PTM639300-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **Aspire Active Air 8 - Alternating Mattress Replacement** is designed for users at high risk of pressure injury, providing advanced pressure care and comfort in both acute and long-term care environments. Equipped with AutoDetect, AutoProtect, and Cell-in-Cell technology, it delivers intelligent automatic adjustment, exceptional stability, and the lowest interface pressures to enhance patient care and safety.

Features

- **Digital Pump:** Easy-to-read LCD screen displays clear cycle times, operation modes, and fault codes for efficient setup and monitoring.
- **AutoLock:** Automatically locks the control panel after 30 seconds to prevent unintentional setting changes.
- **AutoDetect:** Adjusts pressure output automatically based on the user's weight and size for personalised care and protection.
- **AutoProtect:** Increases air pressure automatically when the bed head is raised 45° or more to prevent bottoming out.
- **Continuous Low Pressure (CLP):** Enhances immersion and envelopment for optimal pressure

redistribution and comfort.

- **24-Hour Transport Mode:** Maintains inflation during power outages, ensuring ongoing support and safety.
- **Sleep Quality:** Quiet and smooth pump technology minimises noise and vibration for improved rest.
- **Infection Control:** Welded seams and zip protection reduce the risk of fluid and bacterial ingress.
- **Skin Microclimate:** Multi-stretch cover minimises shear and features a high moisture vapour transfer rate (MVTR) for breathability.
- **Workplace Health & Safety:** Integrated power cord storage eliminates tripping hazards in care environments.
- **Pressure Care Outcomes:** Cell-in-Cell design maximises pressure reduction throughout the alternation cycle.
- **Laundry & Hygiene:** Machine washable cover at 95°C ensures superior cleaning and infection control.

Specifications

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Pressure Care Rating
PTM639300	Active Air 8 Alternating Mattress Replacement – 2000 mm Single	2000 mm	900 mm	200 mm	30 – 200 kg	High Risk
PTM639310	Active Air 8 Alternating Mattress Replacement – 2000 mm King Single	2000 mm	1050 mm	200 mm	30 – 200 kg	High Risk