




## ASPIRE ALTITUDE SINGLE ACTION RECLINER



**SKU:** CHP227705

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Aspire Altitude Recliner is an innovative electric recliner chair crafted to ensure long-lasting comfort. Whether it's for aged care facilities, homes, or hospitals, this recliner seamlessly combines mobility, durability, and user-friendly features.

### FEATURES

- **All-Day Battery Power:** Keep the recliner operational without constant power plug-ins, thanks to its built-in battery. Please note that recharging is necessary.
- **Smooth Relocation:** With its 75 mm Tente castors, easily move the recliner between locations or residents. The rear locking feature enhances safety during movement.
- **Pressure Injury Prevention:** The chair's seating surface is made from a multi-stretch vapour-permeable material, specifically designed to prevent pressure injuries and promote skin health.
- **Medical-Grade Vinyl:** Ensuring extended durability and better infection control, the chair's outer surface is made from medical-grade vinyl.
- **Easy Transfers:** Facilitate safe lateral patient transfers with the chair's removable armrests.
- **Triple-Tiered Backrest:** Enhance the user's comfort with a three-tiered backrest cushion, promoting extended relaxation.
- **Optional Oncology Armrests:** These armrests (CHP227715) are designed to offer wide forearm support, particularly aiding in procedures like blood transfusions. This feature elevates patient comfort and ensures effective treatment administration.
- **Versatile Armrest Table:** The optional armrest table (CHP227720) is a convenient addition for placing personal items or medical instruments during treatments.

The Aspire Altitude Recliner is more than just a chair – it's a comprehensive solution for varied care

needs. It provides an optimal blend of comfort, mobility, and durability, making it an indispensable asset for various health and care facilities. This recliner is tailored to enhance user experience, ensuring safety, comfort, and efficiency in every use-case scenario.

## SPECIFICATIONS

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP227705	Aspire									
	ALTITUDE									
	Single Action Recliner - Black/Blue	800 mm	710 mm	1200 mm	150 kg	520 mm	510 mm	560 mm	560 mm	710 mm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)