



ASPIRE ALUMINIUM FOLDING WALKING FRAME

~~\$107.00~~ Original price was: \$107.00. \$50.00 Current price is: \$50.00.



SKU: WAF706700

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

This lightweight and sturdy frame promises to provide robust support for those needing assistance with ambulation. The frame is made from high-quality aluminium, ensuring it's both durable and easy to manoeuvre. The side-folding feature is a standout, allowing for effortless storage and portability, making it an ideal choice for those on-the-go or with limited storage space. The ergonomic hand grips offer superior comfort, and height adjustment options make this frame highly versatile, accommodating a wide range of user heights.

FEATURES

- Provides assistance for those who have difficulty walking
- Folds flat for easy transportation and storage
- Lightweight aluminium frame ideal for indoor use
- Height adjustable to suit user requirements
- Supplied with four rubber tips
- Width: 610 mm

SPECIFICATIONS

Code	Name	Total Height	Safe Working Load	Unit Weight
WAF706700	Aluminium Folding Walking Frame - Small	710 - 787 mm	127 kg	2 kg