



ASPIRE ALUMINIUM FOLDING WALKING FRAME



SKU: WAF706700

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

This lightweight and sturdy frame promises to provide robust support for those needing assistance with ambulation. The frame is made from high-quality aluminium, ensuring it's both durable and easy to manoeuvre. The side-folding feature is a standout, allowing for effortless storage and portability, making it an ideal choice for those on-the-go or with limited storage space. The ergonomic hand grips offer superior comfort, and height adjustment options make this frame highly versatile, accommodating a wide range of user heights.

FEATURES

- Provides assistance for those who have difficulty walking
- Folds flat for easy transportation and storage
- Lightweight aluminium frame ideal for indoor use
- Height adjustable to suit user requirements
- Supplied with four rubber tips
- Width: 610 mm

SPECIFICATIONS

Code	Name	Total Height	Safe Working Load	Unit Weight
WAF706700	Aluminium Folding Walking Frame - Small	710 - 787 mm	127 kg	2 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au