

ASPIRE ASSIST 3 MANUAL WHEELCHAIR



SKU: MWS452220-H



PRODUCT DESCRIPTION

Purpose & Benefits of the Aspire Assist 3 Manual Wheelchair

The Aspire Assist 3 Manual Wheelchair is a highly versatile and reliable mobility solution designed to meet the needs of individuals requiring lightweight and compact transport. With its adjustable features and durable construction, the Assist 3 ensures optimal comfort, safety, and convenience for both users and caregivers. This wheelchair is ideal for everyday use in home, clinical, or outdoor environments, offering independence and ease of use.

Features

- **Lightweight Aluminium Frame:** Made from durable yet lightweight materials, the wheelchair is easy to manoeuvre and transport without compromising strength.
- **Swing-Away Leg Rests:** The adjustable leg rests can swing away or be removed, making transfers in and out of the wheelchair more convenient.
- **Foldable Design:** The compact folding mechanism allows for easy storage and transport, making it an excellent choice for travel or tight spaces.
- **Height Adjustable Armrests:** Armrests can be adjusted to suit individual user needs, providing enhanced comfort and support.
- **Puncture-Proof Tyres:** Equipped with solid, puncture-proof tyres, ensuring smooth navigation over a variety of surfaces without the need for maintenance.
- **User-Friendly Braking System:** Features an easy-to-operate brake system that enhances safety during use, especially on inclines or uneven terrain.

Specifications

Model	Material	Seat Width	Overall Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)	Armrest Height (mm)	Min. Chair Weight (kg)	SWL (kg)
MWS452200	Dacron	35cm (13.7")	570	400	460 - 500	440	210	17.3	140
MWS452210	Dacron	40cm (15.7")	620	400	460 - 500	440	210	17.5	140
MWS452220	Dacron	45cm (17.7")	670	400	460 - 500	440	210	17.7	140
MWS452230	Dacron	50cm (19.6")	720	400	460 - 500	440	210	17.9	140
MWS452240	WARD Vinyl	35cm (13.7")	570	400	460 - 500	440	210	17.4	140
MWS452250	WARD Vinyl	40cm (15.7")	620	400	460 - 500	440	210	17.6	140
MWS452260	WARD Vinyl	45cm (17.7")	670	400	460 - 500	440	210	17.8	140
MWS452270	WARD Vinyl	50cm (19.6")	720	400	460 - 500	440	210	18	140
MWS452280	HD WARD Vinyl	57cm (22.4")	790	400	460 - 500	440	210	18.9	180
MWS452290	HD WARD Vinyl	61cm (24")	830	400	460 - 500	440	210	19.2	180

The Aspire Assist 3 Wheelchair combines practical features with robust construction, offering users and caregivers a dependable and comfortable mobility solution for everyday needs.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au