

ASPIRE ASSIST 3 MANUAL WHEELCHAIR

\$521.10 - \$822.00 Price range: \$521.10 through \$822.00



SKU: MWS452220-H



PRODUCT DESCRIPTION

Purpose & Benefits of the Aspire Assist 3 Manual Wheelchair

The Aspire Assist 3 Manual Wheelchair is a highly versatile and reliable mobility solution designed to meet the needs of individuals requiring lightweight and compact transport. With its adjustable features and durable construction, the Assist 3 ensures optimal comfort, safety, and convenience for both users and caregivers. This wheelchair is ideal for everyday use in home, clinical, or outdoor environments, offering independence and ease of use.

Features

- **Lightweight Aluminium Frame:** Made from durable yet lightweight materials, the wheelchair is easy to manoeuvre and transport without compromising strength.
- **Swing-Away Leg Rests:** The adjustable leg rests can swing away or be removed, making transfers in and out of the wheelchair more convenient.
- **Foldable Design:** The compact folding mechanism allows for easy storage and transport, making it an excellent choice for travel or tight spaces.
- **Height Adjustable Armrests:** Armrests can be adjusted to suit individual user needs, providing enhanced comfort and support.
- **Puncture-Proof Tyres:** Equipped with solid, puncture-proof tyres, ensuring smooth navigation over a variety of surfaces without the need for maintenance.
- **User-Friendly Braking System:** Features an easy-to-operate brake system that enhances safety during use, especially on inclines or uneven terrain.

Specifications

Model	Material	Seat Width	Overall Width (mm)	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)	Armrest Height (mm)	Min. Chair Weight (kg)	SWL (kg)
MWS452200	Dacron	35cm (13.7")	570	400	460 - 500	440	210	17.3	140
MWS452210	Dacron	40cm (15.7")	620	400	460 - 500	440	210	17.5	140
MWS452220	Dacron	45cm (17.7")	670	400	460 - 500	440	210	17.7	140
MWS452230	Dacron	50cm (19.6")	720	400	460 - 500	440	210	17.9	140
MWS452240	WARD Vinyl	35cm (13.7")	570	400	460 - 500	440	210	17.4	140
MWS452250	WARD Vinyl	40cm (15.7")	620	400	460 - 500	440	210	17.6	140
MWS452260	WARD Vinyl	45cm (17.7")	670	400	460 - 500	440	210	17.8	140
MWS452270	WARD Vinyl	50cm (19.6")	720	400	460 - 500	440	210	18	140
MWS452280	HD WARD Vinyl	57cm (22.4")	790	400	460 - 500	440	210	18.9	180
MWS452290	HD WARD Vinyl	61cm (24")	830	400	460 - 500	440	210	19.2	180

The Aspire Assist 3 Wheelchair combines practical features with robust construction, offering users and caregivers a dependable and comfortable mobility solution for everyday needs.