

ASPIRE BARIATRIC SHOWER COMMODE



SKU: BTC061200

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aspire Bariatric Shower Commode is a specialised bathroom aid designed to provide a safe and convenient toileting solution for individuals with limited mobility. It features a durable frame made of high-strength steel and a weight capacity of up to 350 kg, making it suitable for bariatric users. It also comes with removable swing-away footrests, padded armrests, and a padded backrest for added comfort. The commode bucket is easy to remove and clean, while the locking rear casters ensure stability and safety during transfers.

FEATURES

- Super Heavy Duty Bariatric design with 790 mm seat and 350 kg SWL.
- Anti-bacterial frame coating for improved infection control.
- Removable arms to assist with positioning and allow lateral slide transfers.

SPECIFICATIONS

Code	Name	Seat Width	Safe Working Load
BTC061200	Bariatric Shower Commode	790 mm	350 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au