



ASPIRE BARIATRIC SHOWER COMMODE

\$4,404.00



SKU: BTC061200

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aspire Bariatric Shower Commode is a specialised bathroom aid designed to provide a safe and convenient toileting solution for individuals with limited mobility. It features a durable frame made of high-strength steel and a weight capacity of up to 350 kg, making it suitable for bariatric users. It also comes with removable swing-away footrests, padded armrests, and a padded backrest for added comfort. The commode bucket is easy to remove and clean, while the locking rear casters ensure stability and safety during transfers.

FEATURES

- Super Heavy Duty Bariatric design with 790 mm seat and 350 kg SWL.
- Anti-bacterial frame coating for improved infection control.
- Removable arms to assist with positioning and allow lateral slide transfers.

SPECIFICATIONS

| Code | Name | Seat Width | Safe Working Load |
|-----------|--------------------------|------------|-------------------|
| BTC061200 | Bariatric Shower Commode | 790 mm | 350 kg |