



# **ASPIRE BARIATRIC SHOWER COMMODE**

\$4,026.00 Original price was:

\$4,026.00.<u>\$3,623.00</u>Current price is: \$3,623.00.



**SKU: BTC061200** 

# **PRODUCT DESCRIPTION**

#### **PURPOSE & BENEFITS**

The Aspire Bariatric Shower Commode is a specialised bathroom aid designed to provide a safe and convenient toileting solution for individuals with limited mobility. It features a durable frame made of high-strength steel and a weight capacity of up to 350 kg, making it suitable for bariatric users. It also comes with removable swing-away footrests, padded armrests, and a padded backrest for added comfort. The commode bucket is easy to remove and clean, while the locking rear casters ensure stability and safety during transfers.

### **FEATURES**

- Super Heavy Duty Bariatric design with 790 mm seat and 350 kg SWL.
- Anti-bacterial frame coating for improved infection control.
- Removable arms to assist with positioning and allow lateral slide transfers.

# **SPECIFICATIONS**

Code	Name	<b>Seat Width</b>	Safe Working Load
BTC061200	Bariatric Shower Commode	790 mm	350 kg