




ASPIRE BATH TRANSFER BENCH WITH ROTATING SEAT

\$579.00



SKU: BTS150310

PRODUCT DESCRIPTION

PURPOSE & BENEFITS OF THE ASPIRE BATH TRANSFER BENCH WITH ROTATING SEAT

The Aspire Bath Transfer Bench with Rotating Seat is designed to enhance safety and ease during bathing for individuals with limited mobility. Its rotating seat facilitates smooth transfers into and out of the bathtub, reducing the risk of slips and falls. This bench is ideal for home use, providing users with greater independence and confidence during bathing routines.

Features

- **360° Rotating Seat:** The seat rotates fully, with lockable positions every 90°, accommodating both left and right bath orientations and ensuring secure transfers.
- **Height Adjustable Legs:** Legs adjust in 12.5 mm increments to suit various bath-to-floor heights, offering a customised fit for user comfort.
- **Comfortable Padded Surfaces:** Padded backrest, seat, and arms provide a warmer, more comfortable sitting surface during use.
- **Slip-Resistant Suction Tips:** Equipped with suction rubber tips to enhance stability and prevent movement during transfers.
- **Lightweight Aluminium Construction:** The bench is lightweight for ease of use and portability, without compromising on strength.
- **Safety Belt Included:** Features a safety belt to secure the user during transfers, adding an extra layer of safety.

Specifications

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
BTS150310	Aspire Transfer Bench with Rotating Seat	475 mm	950 mm	795 - 905 mm	136 kg	Aluminium Frame	8.1 kg	345 mm	425 mm	505 - 615 mm	505 - 615 mm

The Aspire Bath Transfer Bench with Rotating Seat offers a practical solution for safer and more comfortable bathing experiences, promoting independence and reducing caregiver strain.