




ASPIRE CLASSIC 6" SEAT WALKER



SKU: WAF709600-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Unveiling the epitome of mobility, comfort, and style – the Aspire Classic 6" Seat Walker. Crafted meticulously with a robust aluminium frame, this walker brings together functionality and user comfort. Whether you're navigating indoor spaces or taking it outdoors, its features promise stability, ease, and a touch of luxury for every step you take.

FEATURES

- **Durable Construction:** Built with a robust aluminium frame that promises longevity while its easy wipe-clean surfaces ensure the walker remains pristine and hygienic.
- **Ultra-soft Hand Brakes:** Thoughtfully designed for users with ailments like arthritis or wrist/hand joint pain, these brakes lock effortlessly. The palm ball design adds an extra layer of comfort during operation.
- **Ergonomically Designed Handles:** Anatomically tailored and latex-free, these handles promise a comfortable grip, making them a boon, especially for those battling arthritis.
- **Personalised Comfort:** Adjustable handle height ensures every user finds their perfect posture, facilitating a comfortable walking experience.
- **Compact Storage:** Its easy folding mechanism ensures the walker can be stowed away conveniently, making transportation and storage a breeze.
- **Secret Storage:** The padded seat not only offers a comfortable resting place but also flips up to reveal a discreet pocket for essentials. Additionally, a storage bag is provided to keep personal items secure.
- **Seamless Mobility:** Equipped with solid, puncture-proof 6" wheels, the walker offers smooth manoeuvrability, especially indoors, ensuring users can navigate spaces with ease.

SPECIFICATIONS

Overall Height	760 - 915 mm
Overall Width	615 mm
Overall Length	625 mm
Seat Depth	305 mm
Seat Width	355 mm
Seat Height	525 mm
Backrest Height	770 mm
Folded Width	200 mm
Wheel Size	6"
Unit Weight	6 kg
Safe Working Load	130 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au