



ASPIRE CLASSIC 8" SEAT WALKER

\$201.00



SKU: WAF705700-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Discover an unparalleled mobility experience with the Aspire Classic 8" Seat Walker. Seamlessly combining ergonomic design, utility, and aesthetics, this walker emerges as an essential companion for those seeking enhanced mobility. Made with a resilient aluminium frame and tailored to cater to individual comfort, it's clear that the walker is designed to keep the user at its core.

FEATURES

- **Sturdy Construction:** Crafted with a robust aluminium frame, the walker assures long-term durability. The easy wipe-clean surfaces make maintenance both convenient and swift, ensuring your walker remains in pristine condition.
- **Ultra-soft Hand Brakes:** Designed with the user's comfort in mind, the brakes, featuring a palm ball, lock in seamlessly. They're particularly beneficial for those with conditions like arthritis or wrist/hand joint pain, offering effortless operation without compromising on security.
- **Ergonomic Handles:** Shaped anatomically and devoid of latex, these handles are built for an unparalleled grip. They cater especially to those with arthritis, promising comfort with every use.
- **Customized Comfort:** The walker's adjustable handle height allows users to find their optimal posture, ensuring comfort during use.
- **Compact Design:** A simple folding mechanism ensures that storing and transporting the walker is hassle-free, making it a perfect travel companion.
- **Integrated Storage:** A padded seat, not only offers a comfortable rest spot but also flips up to unveil a discreet pocket for essentials. An additional storage bag ensures all your belongings remain secure.

- **Enhanced Mobility:** With solid, puncture-proof 8" wheels, this walker is a mobility dream, especially outdoors. The larger wheel size ensures smooth manoeuvrability on varying terrains.

The Aspire Classic 8" Seat Walker stands out as a beacon of modern mobility. Every feature, from its sturdy frame to its user-centric design, screams of sophistication and thoughtfulness. Whether you're navigating the corridors of your home or the pathways of a park, this walker ensures a smooth, safe, and comfortable journey. Step out in confidence and style with this exceptional mobility aid that promises to transform everyday walks into delightful journeys.

SPECIFICATIONS

Code	Display Name	Seat Width	Total Height	Seat To Floor Height	Safe Working Load	Unit Weight	Colour
WAF705700RE	Aspire Classic 8 Seat Walker - Red	355 mm	817 - 972 mm	580 mm	130 kg	7 kg	Red
WAF705700BL	Aspire Classic 8 Seat Walker - Blue	355 mm	817 - 972 mm	580 mm	130 kg	7 kg	Blue