



ASPIRE CONFIMOTION ACTIVCARE BED

\$3,747.00 - \$5,952.00 Price range: \$3,747.00 through \$5,952.00



SKU: BEB046740-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **Aspire ComfiMotion ActivCare Bed** combines advanced functionality, clinical positioning, and contemporary design to provide comfort, safety, and independence for home and aged care users. Designed for both user wellbeing and caregiver efficiency, it offers a range of adjustable features to support mobility, circulation, and relaxation while maintaining a stylish appearance suitable for any home environment.

Features

- **Ultra Low Height:** Platform lowers to 195 mm to reduce fall risk and enhance user safety.
- **Customisable Entry/Exit Height:** Programmable single-button operation enables personalised height settings for safe transfers.
- **Under Bed Clearance:** Allows easy access for patient lifters and overbed tables.
- **Battery Backup:** Ensures continued function and safety during power outages.
- **Night Assist:** Luminous underbed halo light provides visibility for safe movement at night.
- **Illuminated Handset:** Large, easy-to-read interface with intuitive controls.

- **Massage Function:** Adjustable vibration modes offer neck, shoulder, and back relief.
- **Hi-Lo Functionality:** Simplifies manual handling and supports safer patient care.
- **Back Raise & Head Tilt:** Promotes postural support and comfort for reading, TV, or relief of respiratory and reflux conditions.
- **Auto-Regression Backrest:** Prevents user sliding, improving comfort and pressure redistribution.
- **Clinical Positioning:** One-touch function raises back and legs to reduce cardiac strain and improve rest.
- **Leg Raise & Zero Gravity:** Enhances circulation and relieves lower back tension and fatigue.
- **Memory Function:** Saves preferred positions for one-touch recall.
- **Seated Position:** Assists with upright dining and bed transfers.
- **Trendelenburg & Reverse:** Supports repositioning under professional supervision.
- **Lockable Castors:** Eight heavy-duty dual castors ensure secure stability and easy movement.

Specifications

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Colour
BEB046740	Aspire ComfiMotion ActivCare Bed – Long Single (Headboard Separate)	2030 mm	920 mm	195 – 700 mm	250 kg	Portobello Grey
BEB046755	Aspire ComfiMotion ActivCare Bed – King Single (Headboard Separate)	2030 mm	1070 mm	200 – 775 mm	250 kg	Portobello Grey

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Colour
BEB046760	Aspire ComfiMotion ActivCare Bed – Long Double (Headboard Separate)	2030 mm	1370 mm	200 – 675 mm	360 kg	Portobello Grey
BEB046770	Aspire ComfiMotion ActivCare Bed – Queen (Headboard Separate)	2030 mm	1520 mm	200 – 675 mm	360 kg	Portobello Grey