

ASPIRE CONFIMOTION BED COLLECTION - BED RAIL



SKU: BEA830210

PRODUCT DESCRIPTION

Purpose & Benefits

The **Aspire ComfiMotion Bed Collection - Bed Rail** enhances user safety and comfort by improving bed edge awareness and providing multiple grip points for repositioning and mobility. Designed for compatibility with the ComfiMotion bed range, it assists users in safely entering, exiting, and adjusting positions while in bed.

Features

- Improves bed edge awareness for enhanced user safety
- Provides multiple grip points to assist with repositioning and mobility
- Compatible with the Aspire ComfiMotion Bed Collection
- Extendable version (BEA830205) allows the rail to expand for fall prevention and contract for easier bed access
- Durable construction for stability and long-term reliability

Specifications

- **Model Code:** BEA830205 (Extendable Version)

- **Height from Mattress Deck:** 565 mm
- **Width of Rail:** 685 mm
- **Compatibility:** Aspire ComfiMotion Bed Collection

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au