



ASPIRE DELUXE BEDSIDE COMMODORE



SKU: BEC04511-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aspire Deluxe Bedside Commode is a thoughtful integration of comfort and functionality, designed for users who find it challenging to move between the bedroom and bathroom. With its easy-to-clean surfaces and adjustable features, it ensures both convenience and hygiene, making daily life a bit easier for those in need.

FEATURES

- **Bedroom Toileting Solution:** Specifically designed for those who can't frequently travel between the bedroom and bathroom, ensuring they have a comfortable and accessible toileting solution close at hand.
- **Colour Options:** Available in two aesthetically pleasing colours - Champagne or Slate Vinyl, allowing users to choose one that best suits their personal preference or bedroom decor.
- **Easy Clean Vinyl Surfaces:** The use of vinyl surfaces ensures that the commode is easy to clean, maintaining hygiene with minimal effort.
- **Removable Bowl with Added Features:** The bowl can be effortlessly removed for cleaning, and it comes with a lid for hygiene and a carry handle for easy transportation.
- **Customized Comfort:** The commode's height is adjustable, ensuring it can be tailored to fit the user's requirements, offering optimal comfort.
- **Ergonomic Design:** Raised, padded armrests provide additional support, making getting on and off the commode easier and safer. The fold-away, padded seat ensures both comfort and space-saving convenience.
- **Stability and Safety:** Equipped with non-slip tips, the commode guarantees stability, ensuring

users feel safe and secure during use.

The Aspire Deluxe Bedside Commode is more than just a convenience; it's a testament to understanding the needs of those who require extra support in their daily routines.

SPECIFICATIONS

Code	Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
BEC045115	Aspire Deluxe Bedside Commode - Champagne	450 mm	460 mm	410 - 565 mm	160 kg	12 kg
BEC045116	Aspire Deluxe Bedside Commode - Slate	450 mm	460 mm	410 - 565 mm	160 kg	12 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au