



ASPIRE EVOKE 2 WHEELCHAIR



SKU: MWS4495-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aspire Evoke 2 Wheelchair is available in 4 different seat widths. This is a manual wheelchair that is highly adjustable and lightweight for easy transportation. The quick release 24 inch rear wheels will reduce the weight of the chair to assist with even easier transportation. It has various features that make it a stand out option for a wheelchair.

The swing in/out footrests have a design that ensures ultimate comfort ability. Being height adjustable and removable, it allows the user to customise according to their personal requirements.

FEATURES

- The Evoke 2 and Evoke 2 HD are designed for the user who wants a lightweight, highly adjustable manual wheelchair.
- Swing away and removable arm rests this facilitates user transfers and also makes for easier transport and storage.
- Swing-in / out, height adjustable and removable foot rests this facilitates user transfers and allows optimal positioning for the changing needs.
- Strong, Nylon upholstery with a comfortable cushion and padded backrest. For increased comfort, a padded overlay is available as an accessory.
- Depth adjustable seating with rear angle adjustability for more accurate positioning, customisability and increased user comfort.
- Quick release 24" solid rear tyres the wheels are easily removed to reduce the weight of the chair for transport and storage.



SPECIFICATIONS

Model	14"	16"	18"	20"	22"	24"
Seat Width:	35cm	40cm	45cm	50cm	57cm	61cm
Seat Depth:	40 to 47cm	s40 to 47cm	s 40 to 47cr	m40 to 47cr	m40 to 47cm	n 40 to 47cm
Back Height:	43 to 49cm	43 to 49cm	43 to 49cr	n43 to 49cr	m43 to 49cm	n 43 to 49cm
Seat to Floor Height:	44 to 48cm	s44 to 48cm	s 44 to 48cr	n44 to 48cr	n 44 to 48cr	m44 to 48cm
Weight:	14.8kg	14.9kg	15kg	15.1kg	16kg	16.2kg
Weight Capacity:	140kg	140kg	140kg	140kg	180kg	180kg
Warranty:	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years

LIGHTWEIGHT ASPIRE WHEELCHAIR

The Evoke 2 is intended for a person that requires an adjustable lightweight manual wheelchair. The Contoured Padded Overlay is a stylish favoured accessory for the Aspire Evoke 2.

Buy the Aspire Evoke Wheelchair Australia online at MediMart using Afterpay or NDIS method today!

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom
343 Brighton Rd, North
Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777 elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771 reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT



Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au