




## **ASPIRE FLEX ADJUSTABLE SEAT WALKER 6" WHEEL**

**\$205.00**



**SKU: WAF709660CH**

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

Enhance your mobility with the Aspire Flex Adjustable Seat Walker, featuring 6" wheels designed for smooth navigation and stability. This versatile mobility aid combines the convenience of a walker with the comfort of a seat, allowing users to rest as needed during their activities.

### **FEATURES**

- **Adjustable Height:** The handles of the walker are adjustable, ensuring a comfortable fit for users of various heights.
- **Comfortable Seating:** A padded seat provides a convenient and comfortable place to rest when not in motion.
- **Compact Design:** The walker is engineered to be easily folded, making it simple to store or transport.
- **Durable Construction:** Built with high-quality materials, it's made to withstand daily use and provide reliable support.
- **Safety Brakes:** Equipped with easy-to-use brakes for added safety and control when walking or sitting.
- **Lightweight Frame:** The frame is constructed to be light yet sturdy, enhancing maneuverability without sacrificing strength.

Whether navigating indoor spaces, traversing outdoor terrain, or needing a quick rest, the Aspire Flex Adjustable Seat Walker with 6" wheels is an excellent choice for individuals seeking a reliable walking aid that doesn't compromise on comfort or convenience.

## SPECIFICATIONS

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Width	Seat To Floor Height	Seat To Footrest
WAF709660CH	Aspire									
	Flex									
	Adjustable Seat	730 mm	570 mm	730 - 830 mm	136 kg	6.3 kg	Champagne	360 mm	430 - 520 mm	430 - 520 mm
	Walker 6" Wheel									