



## ASPIRE HIGH BACK CLASSIC DAY CHAIR



**SKU:** CHP20805-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Elevate comfort and support to new heights with the Aspire High Back Classic Day Chair. Designed with orthopaedic considerations in mind, this chair melds functionality with unparalleled support, ensuring users enjoy comfort throughout the day. Its high-back contouring, plush padding, and adjustability make it a must-have for those seeking a combination of style, durability, and comfort.

### FEATURES

- **Contoured Backrest:** The high back design isn't just about aesthetics—it's meticulously crafted to offer increased postural support, ensuring users maintain a healthy posture even after prolonged seating.
- **Luxurious Comfort:** From the softly padded armrests to the seat and back upholstery, every inch of this chair spells comfort. It's the perfect blend of firm support and soft cushioning.
- **Hassle-Free Maintenance:** The easy clean vinyl surfaces ensure that tidying up is a breeze, making it perfect for everyday use.
- **Customizable Height:** With an easy push-button adjustment feature, users can set the chair to their preferred height, making it adaptable to changing needs.
- **Sturdy yet Light:** Despite its generous weight capacity of 160 kg, the chair remains lightweight, ensuring it can be moved around with ease.
- **Durable Design:** The powder-coated steel frame isn't just robust; it's also easy to clean, ensuring the chair remains in pristine condition for years.
- **Safety First:** The non-slip tips on the chair legs ensure utmost stability, providing peace of mind during sitting or transferring sessions.
- **Colour Variants:** Depending on individual preferences, users can opt for either the

sophisticated Slate or the elegant Champagne Vinyl finish.

The Aspire High Back Classic Day Chair goes beyond being a mere seating solution—it's a testament to ergonomic design, comfort, and adaptability. Whether you're settling in for a leisurely read or looking for orthopaedic support, this chair ensures that every moment spent on it is one of absolute comfort. Its durable design, combined with its aesthetic appeal, makes it a worthy addition to any living space.

## SPECIFICATIONS

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Safe Working Load	Colour
CHP208050	Aspire High Back Classic Day Chair - Champagne Vinyl	490 mm	490 mm	720 mm	425 - 610 mm	160 kg	Champagne Vinyl
CHP208055	Aspire High Back Classic Day Chair - Slate Vinyl	490 mm	490 mm	720 mm	425 - 610 mm	160 kg	Slate Vinyl

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)