


## ASPIRE HOMECARE OVER TOILET AID



**SKU:** BTT145025

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The **Aspire Homecare Over Toilet Aid** is designed to enhance safety, comfort, and independence during toileting. Featuring adjustable height and sturdy armrests, it assists users in sitting and standing with greater ease. Ideal for individuals with reduced mobility or balance challenges, this aid promotes confidence and stability while reducing the risk of falls in the bathroom.

### FEATURES

- Adjustable height to accommodate a wide range of users
- Durable, non-slip high-grip feet ensure stability and safety
- Raised armrests provide additional support for safe transfers
- Lightweight yet sturdy frame for everyday use
- Easy to clean and maintain for long-lasting hygiene

### SPECIFICATIONS

- **Height Range:** 595 mm – 750 mm
- **Total Width:** 520 mm
- **Total Depth:** 480 mm
- **Seat Depth:** 410 mm
- **Seat Width:** 340 mm
- **Seat to Floor Height:** 405 mm – 560 mm
- **Width Between Arms:** 475 mm
- **Safe Working Load (SWL):** 100 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)