



## **ASPIRE HOMECARE OVER TOILET AID**

\$128.00 Original price was: \$128.00.\$115.20 Current

price is: \$115.20.



SKU: BTT145025

## PRODUCT DESCRIPTION

#### **PURPOSE & BENEFITS**

The **Aspire Homecare Over Toilet Aid** is designed to enhance safety, comfort, and independence during toileting. Featuring adjustable height and sturdy armrests, it assists users in sitting and standing with greater ease. Ideal for individuals with reduced mobility or balance challenges, this aid promotes confidence and stability while reducing the risk of falls in the bathroom.

### **FEATURES**

- Adjustable height to accommodate a wide range of users
- Durable, non-slip high-grip feet ensure stability and safety
- Raised armrests provide additional support for safe transfers
- Lightweight yet sturdy frame for everyday use
- Easy to clean and maintain for long-lasting hygiene

# **SPECIFICATIONS**

• Height Range: 595 mm - 750 mm

• Total Width: 520 mm • Total Depth: 480 mm • Seat Depth: 410 mm • Seat Width: 340 mm

• Seat to Floor Height: 405 mm - 560 mm

• Width Between Arms: 475 mm • Safe Working Load (SWL): 100 kg