



ASPIRE KITCHEN/PERCHING STOOL

\$256.00 Original price was: \$256.00.<u>\$230.40</u>Current price is: \$230.40.



SKU: N/A

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Aspire Kitchen/Perching Stool** is designed to provide comfort, safety, and support for users who need to rest while performing kitchen or household tasks. Featuring an anterior tilt to assist with standing and sitting, this stool promotes independence and reduces strain on the legs and lower back. Ideal for elderly users, individuals with limited mobility, or those recovering from injury, it offers both practicality and comfort in everyday use.

FEATURES

- Anterior tilt seat aids the sitting-to-standing process for improved independence
- Height adjustable with easy push-button system to suit individual needs
- Padded seat and backrest provide excellent comfort and support
- Non-slip rubber tips enhance stability and user safety
- Powder-coated steel frame offers strength, durability, and easy cleaning
- Compact design ideal for kitchen or home use

SPECIFICATIONS

Model Code: CHP213500

• Display Name: Aspire Kitchen/Perching Stool - Champagne

Seat Depth: 310 mmSeat Width: 375 mm

• Backrest Height: 400 mm

Seat to Floor Height: 540 - 690 mm
Safe Working Load (SWL): 175 kg



• Frame Colour: Champagne Vinyl