




## ASPIRE KITCHEN/PERCHING STOOL



**SKU:** CHP213500

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The **Aspire Kitchen/Perching Stool** is designed to provide comfort, safety, and support for users who need to rest while performing kitchen or household tasks. Featuring an anterior tilt to assist with standing and sitting, this stool promotes independence and reduces strain on the legs and lower back. Ideal for elderly users, individuals with limited mobility, or those recovering from injury, it offers both practicality and comfort in everyday use.

### FEATURES

- Anterior tilt seat aids the sitting-to-standing process for improved independence
- Height adjustable with easy push-button system to suit individual needs
- Padded seat and backrest provide excellent comfort and support
- Non-slip rubber tips enhance stability and user safety
- Powder-coated steel frame offers strength, durability, and easy cleaning
- Compact design ideal for kitchen or home use

### SPECIFICATIONS

- **Model Code:** CHP213500
- **Display Name:** Aspire Kitchen/Perching Stool - Champagne
- **Seat Depth:** 310 mm
- **Seat Width:** 375 mm
- **Backrest Height:** 400 mm
- **Seat to Floor Height:** 540 – 690 mm
- **Safe Working Load (SWL):** 175 kg

- **Frame Colour:** Champagne Vinyl

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriiffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)