


ASPIRE LIFESTYLE COMMUNITY BED

\$3,463.00



SKU: BEB046040



PRODUCT DESCRIPTION

Purpose & Benefits of the Aspire Lifestyle Community Bed

The **Aspire Lifestyle Community Bed** is designed to provide exceptional comfort, flexibility, and support for individuals in residential and care settings. Featuring modern aesthetics and advanced functionality, this bed is ideal for enhancing quality of life while catering to a wide range of mobility and health needs. It promotes independence and ensures ease of use for both users and caregivers.

Key Features

- **Height Adjustability:** The bed offers adjustable height settings, allowing for safe and comfortable transfers, as well as easy caregiver access during personal care routines.
- **Customisable Positioning:** With motorised controls for head, foot, and height adjustments, the bed enables users to find their ideal position for sleeping, reading, or resting, improving comfort and circulation.
- **Modern and Aesthetic Design:** Blending functionality with contemporary style, the Aspire Lifestyle Community Bed fits seamlessly into any home or care environment, offering a sleek and modern appearance.
- **Durable and Safe:** Built with high-quality materials and featuring reinforced construction, the bed is designed for long-lasting performance, ensuring stability and safety for a wide range of users.
- **Easy-to-Use Controls:** The intuitive remote control allows for simple adjustments of bed positions, giving users and caregivers full control over comfort settings at the touch of a button.
- **Split-Side Rail Options:** Optional side rails provide added safety and peace of mind, especially for users who are at risk of falls, while allowing easy access when required.

Specifications

| Specification | Single | King Single |
|----------------------------------|---------------------------|----------------------------|
| Total Length | 2220 mm | 2220 mm |
| Total Length Extended | 2370 mm | 2370 mm |
| Total Width | 995 mm | 1155 mm |
| Height | 245 mm - 645 mm | 245 mm - 645 mm |
| Backrest Incline | 70° | 70° |
| Knee Bend | 30° | 30° |
| Recommended Mattress Size | 2000 mm x 900 mm x 150 mm | 2030 mm x 1050 mm x 150 mm |
| Colour | Gunmetal | Gunmetal |
| Maximum Patient Weight | 180 kg | 180 kg |
| SWL | 220 kg | 220 kg |

The **Aspire Lifestyle Community Bed** is an ideal solution for those seeking a balance of comfort, safety, and independence. Its user-friendly design and modern look make it a versatile choice for care settings or home use, providing exceptional support and flexibility for all users.