


ASPIRE LIFESTYLE COMMUNITY BED

\$3,159.00



SKU: BEB046040



PRODUCT DESCRIPTION

Purpose & Benefits of the Aspire Lifestyle Community Bed

The **Aspire Lifestyle Community Bed** is designed to provide exceptional comfort, flexibility, and support for individuals in residential and care settings. Featuring modern aesthetics and advanced functionality, this bed is ideal for enhancing quality of life while catering to a wide range of mobility and health needs. It promotes independence and ensures ease of use for both users and caregivers.

Key Features

- **Height Adjustability:** The bed offers adjustable height settings, allowing for safe and comfortable transfers, as well as easy caregiver access during personal care routines.
- **Customisable Positioning:** With motorised controls for head, foot, and height adjustments, the bed enables users to find their ideal position for sleeping, reading, or resting, improving comfort and circulation.
- **Modern and Aesthetic Design:** Blending functionality with contemporary style, the Aspire Lifestyle Community Bed fits seamlessly into any home or care environment, offering a sleek and modern appearance.
- **Durable and Safe:** Built with high-quality materials and featuring reinforced construction, the bed is designed for long-lasting performance, ensuring stability and safety for a wide range of users.
- **Easy-to-Use Controls:** The intuitive remote control allows for simple adjustments of bed positions, giving users and caregivers full control over comfort settings at the touch of a button.
- **Split-Side Rail Options:** Optional side rails provide added safety and peace of mind, especially for users who are at risk of falls, while allowing easy access when required.

Specifications

Specification	Single	King Single
Total Length	2220 mm	2220 mm
Total Length Extended	2370 mm	2370 mm
Total Width	995 mm	1155 mm
Height	245 mm - 645 mm	245 mm - 645 mm
Backrest Incline	70°	70°
Knee Bend	30°	30°
Recommended Mattress Size	2000 mm x 900 mm x 150 mm	2030 mm x 1050 mm x 150 mm
Colour	Gunmetal	Gunmetal
Maximum Patient Weight	180 kg	180 kg
SWL	220 kg	220 kg

The **Aspire Lifestyle Community Bed** is an ideal solution for those seeking a balance of comfort, safety, and independence. Its user-friendly design and modern look make it a versatile choice for care settings or home use, providing exceptional support and flexibility for all users.