



ASPIRE LIGHTWEIGHT BOOT SCOOTER

\$3,744.00 Original price was:

\$3,744.00.<u>\$3,369.60</u>Current price is: \$3,369.60.



SKU: SCT670616

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Ultimate** Portable Travel Scooter! Every day is a new adventure with the Boot Scooter Supalite, an exceptionally stable and manoeuvrable 4-wheel travel scooter. Benefiting from a class-leading, ultra-lightweight design and a maximum speed up to 6.5 km/h, the Boot Scooter Supalite keeps up with your active life!

FEATURES

- An exceptionally stable and manoeuvrable 4-wheel travel scooter
- Benefiting from a class-leading, ultra-lightweight design and a maximum speed up to 6.5 km/h
- Compact, foldable seat with armrests
- Easy-to-grip tiller adjustment knob and handles
- Stable, front frame-mounted seat post
- A standard front basket provides storage
- Easy to transport at home, hotels, or airports
- Disassembles into 6 lightweight pieces for easy transport and storage! (heaviest piece 10.7kg)
- Up to 8kg lighter than comparable travel scooters

SPECIFICATIONS

Total Length	Total Width	Safe Working Load	Unit Weight	Maximum Speed	Range
1022 mm	490 mm	120 kg	34 kg	6.5 km/h	18.5km (User Weight: 90.7kg)



Buy This Aspire Lightweight Boot Scooter online at MediMart using Afterpay or <u>NDIS method</u> today! Expert <u>mobility equipment suppliers</u> in Australia.