


## ASPIRE LIGHTWEIGHT BOOT SCOOTER



**SKU:** SCT670616



## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The **Ultimate** Portable Travel Scooter! Every day is a new adventure with the Boot Scooter Supalite, an exceptionally stable and manoeuvrable 4-wheel travel scooter. Benefiting from a class-leading, ultra-lightweight design and a maximum speed up to 6.5 km/h, the Boot Scooter Supalite keeps up with your active life!

### FEATURES

- An exceptionally stable and manoeuvrable 4-wheel travel scooter
- Benefiting from a class-leading, ultra-lightweight design and a maximum speed up to 6.5 km/h
- Compact, foldable seat with armrests
- Easy-to-grip tiller adjustment knob and handles
- Stable, front frame-mounted seat post
- A standard front basket provides storage
- Easy to transport at home, hotels, or airports
- Disassembles into 6 lightweight pieces for easy transport and storage! (heaviest piece 10.7kg)
- Up to 8kg lighter than comparable travel scooters

### SPECIFICATIONS

Total Length	Total Width	Safe Working Load	Unit Weight	Maximum Speed	Range
1022 mm	490 mm	120 kg	34 kg	6.5 km/h	18.5km (User Weight: 90.7kg)

Buy This Aspire Lightweight Boot Scooter online at MediMart using Afterpay or [NDIS method](#) today!  
Expert [mobility equipment suppliers](#) in Australia.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)