



ASPIRE LITE TRANSIT WHEELCHAIR

~~\$566.00~~ Original price was: \$566.00. \$557.10 Current price is: \$557.10.



SKU: MWS442301

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Meet the Aspire Lite Transit Wheelchair - a fusion of lightweight construction and practical design, tailored specifically for the occasional user. Featuring attendant brakes and other user-centric features, this wheelchair offers enhanced mobility and comfort for both the user and the attendant. Its design ensures that each trip, whether short or long, is made in utmost ease and safety.

FEATURES

- **Ultra-Lightweight Construction:** The Aspire Lite Transit Wheelchair boasts an ultra-light folding frame, weighing only 11.5 kg. This ensures it's not just easy to transport but also a cinch to store.
- **Attendant-Controlled Brakes:** The wheelchair is equipped with attendant brakes as a standard feature. This allows for improved maneuverability, especially during turns or stops, giving the attendant better control.
- **Ergonomic Design:** Keeping the attendant's comfort in mind, the wheelchair features ergonomic handles. This ensures that pushing the wheelchair remains an effortless task without causing any discomfort.
- **User-Centric Adjustability:** The push-button, height-adjustable armrests have been specially designed considering users with weak or arthritic hands, allowing them to find their perfect comfort level with ease.
- **Ease of Transfers:** To ensure that getting in and out of the wheelchair remains a hassle-free task, it comes with swing-away, flip-up footplates.
- **Enhanced Terrain Adaptability:** Equipped with 12" rear wheels and easily accessible brakes, the wheelchair promises improved performance over outdoor terrains, ensuring rides remain

smooth regardless of the surface.

The Aspire Lite Transit Wheelchair is a perfect blend of innovation and user-centric design. Its lightweight nature, combined with features like attendant brakes and ergonomic handles, makes it a go-to choice for those seeking reliability and comfort in a package that's easy to manage. Whether it's a trip to the park or a day out in the city, this wheelchair ensures that mobility remains as simple and comfortable as possible.

SPECIFICATIONS

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Safe Working Load	Unit Weight	Colour	Material
MWS442301	Aspire Lite Transit	400 mm	460 mm	465 mm	460 mm	120 kg	11.6 kg	Champagne	Aluminium