

ASPIRE LITE WHEELCHAIR



SKU: MWS442300



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aspire Lite Wheelchair stands out as an epitome of convenience and comfort. Specifically tailored for occasional users, its ultra-lightweight frame ensures easy maneuverability without compromising on stability. It brings together advanced ergonomics and ease of use, ensuring both the user and the caregiver experience maximum comfort and convenience.

FEATURES

- **Ultra-Light Frame:** The hallmark of the Aspire Lite Wheelchair is its ultra-light folding frame. Weighing a mere 9.5 kg, this wheelchair not only ensures easy transportation but also hassle-free storage.
- **Caregiver-Friendly Design:** With ergonomic handles tailored for the caregiver, pushing and maneuvering the wheelchair becomes a breeze, eliminating any undue strain.
- **Adjustable Arm Rests:** Keeping user comfort in the spotlight, the wheelchair comes with push-button, height-adjustable armrests. They are particularly suited for individuals with weak or arthritic hands, ensuring they can adjust it according to their comfort.
- **Facilitating Easy Transfers:** With swing-away, flip-up footplates, the wheelchair ensures users can effortlessly transfer in and out of the chair, enhancing user autonomy.
- **Safety Assured:** The rear wheel locks act as a steadfast safety feature, ensuring stability and security, and giving users and caregivers peace of mind.

The Aspire Lite Wheelchair is a testament to how advanced engineering can elevate the comfort and convenience of wheelchair users and their caregivers. Its design, focused on ease of use and user-friendly features, makes it an indispensable asset for those who seek a blend of lightweight convenience and sturdy reliability. Whether you're using it for short trips or occasional outings, this

wheelchair promises to deliver both functionality and comfort.

SPECIFICATIONS

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Safe Working Load	Unit Weight	Colour	Material
MWS442300	Aspire Lite	400 mm	460 mm	465 mm	460 mm	120 kg	9.5 kg	Blue	Aluminium

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au