

## ASPIRE LOW BACK CLASSIC DAY CHAIR

**\$419.00**



**SKU:** CHP20800-H



## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Aspire Low Back Classic Day Chair is the embodiment of ergonomic design blended with versatility. Built specifically for those seeking postural support without compromising on comfort, this orthopaedic chair offers a perfect seating solution. Its array of user-centric features ensures that each sitting experience is comfortable, safe, and tailored to individual needs.

### FEATURES

- **Low Back, Contoured Design:** The chair's contoured backrest is uniquely designed to offer optimum postural support, ensuring users maintain a comfortable and healthy posture.
- **Premium Comfort:** With softly padded armrests, seat, and back upholstery, users can expect a supremely comfortable seating experience every time.
- **Stable & Secure:** Fitted with non-slip rubber tips, the chair provides unmatched stability, ensuring users can sit or transfer with confidence.
- **Easy Maintenance:** The easy-clean vinyl surfaces ensure that keeping the chair spotless is a breeze. Plus, with choices between Slate and Champagne Vinyl, users can pick a style that best matches their aesthetic.
- **Customisable Height:** The chair comes equipped with an easy push-button adjustment mechanism. This allows users to adjust the chair's height to their preference, ensuring maximum comfort.
- **Mobility & Strength:** While the chair's lightweight design ensures easy maneuverability, its powder-coated steel frame guarantees strength and durability, promising long-term use without wear and tear.

The Aspire Low Back Classic Day Chair is more than just a seating solution; it's a commitment to

comfort, stability, and user-focused design. Whether it's the contoured backrest, adjustable height, or the plush padded surfaces, every feature of this chair is crafted keeping the user's needs at its core. Choose this chair for a harmonious blend of comfort, style, and utility.

## SPECIFICATIONS

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Safe Working Load	Colour
CHP208000	Aspire Low Back Classic Day Chair - Champagne Vinyl	440 mm	460 mm	440 mm	430 - 605 mm	160 kg	Champagne Vinyl
	Aspire Low Back Classic Day Chair - Slate Vinyl	440 mm	460 mm	440 mm	430 - 605 mm	160 kg	Slate Vinyl