




## ASPIRE MAXI ADJUSTABLE BEDSIDE COMMODORE

**\$1,008.00**



**SKU:** BEC045120

## PRODUCT DESCRIPTION

### Purpose & Benefits of the Aspire Maxi Adjustable Bedside Commode

The Aspire Maxi Adjustable Bedside Commode is a versatile and supportive solution for individuals who require toileting assistance near their bedside. Designed to promote dignity and independence, this commode is ideal for those with limited mobility or recovering from surgery. Its height-adjustable frame and generous weight capacity ensure comfort and stability for a range of users, while the easy-to-clean surfaces make it suitable for home care and clinical environments.

### Key Features

- **Height Adjustable Frame:** Customise the height to suit individual needs, providing optimal comfort and ease of access for each user.
- **Generous Weight Capacity:** With a robust design, the Aspire Maxi can accommodate users up to 250 kg, offering secure and reliable support.
- **Removable Padded Seat and Lid:** Features a comfortable padded seat that doubles as a lid, allowing the commode to blend discreetly into a room when not in use.
- **Easy-to-Clean Surfaces:** Made from materials that are simple to wipe down, ensuring hygiene is maintained with minimal effort.
- **Non-Slip Rubber Feet:** Equipped with rubber-tipped feet to provide stability on various flooring types, reducing the risk of slipping during use.

### Specifications

- **Seat Width:** Adjustable from 57 cm to 77 cm
- **Height Range:** Adjustable from 45 cm to 57 cm
- **Safe Working Load:** 310 kg

- **Material:** Powder-coated steel frame with padded seat and backrest
- **Cleaning:** Wipe clean with a damp cloth or disinfectant