




## ASPIRE MAXI ADJUSTABLE BEDSIDE COMMODOE



**SKU:** BEC045125

## PRODUCT DESCRIPTION

### Purpose & Benefits of the Aspire Maxi Adjustable Bedside Commode

The Aspire Maxi Adjustable Bedside Commode is a versatile and supportive solution for individuals who require toileting assistance near their bedside. Designed to promote dignity and independence, this commode is ideal for those with limited mobility or recovering from surgery. Its height-adjustable frame and generous weight capacity ensure comfort and stability for a range of users, while the easy-to-clean surfaces make it suitable for home care and clinical environments.

### Key Features

- **Height Adjustable Frame:** Customise the height to suit individual needs, providing optimal comfort and ease of access for each user.
- **Generous Weight Capacity:** With a robust design, the Aspire Maxi can accommodate users up to 250 kg, offering secure and reliable support.
- **Removable Padded Seat and Lid:** Features a comfortable padded seat that doubles as a lid, allowing the commode to blend discreetly into a room when not in use.
- **Easy-to-Clean Surfaces:** Made from materials that are simple to wipe down, ensuring hygiene is maintained with minimal effort.
- **Non-Slip Rubber Feet:** Equipped with rubber-tipped feet to provide stability on various flooring types, reducing the risk of slipping during use.

### Specifications

- **Seat Width:** Adjustable from 57 cm to 77 cm
- **Height Range:** Adjustable from 45 cm to 57 cm
- **Safe Working Load:** 310 kg

- **Material:** Powder-coated steel frame with padded seat and backrest
- **Cleaning:** Wipe clean with a damp cloth or disinfectant

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)