



ASPIRE MAXI DAY CHAIR - UTILITY CHAIR

\$1,277.00



SKU: CHP208010

PRODUCT DESCRIPTION

This sturdy bariatric seating support surface is an ideal bariatric seating solution with adjustable seat width, depth & backrest height.

FEATURES

- Generous 310kg safe working load, a true bariatric solution.
- Ergonomic seat with anterior tilt to improve posture and reduce fatigue.
- Height, width and depth adjustability in seat and back for customised comfort.
- Easily manoeuvrable with push handle and transport wheels.
- Easy clean, breathable and comfortable vinyl upholstery and powder coated frame for effortless maintenance and longevity.

SPECIFICATIONS

Name	Total Height	Safe Working Load	Material	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
Maxi Day Chair	1120 - 1440 mm	310 kg	Vinyl	45 kg	390 - 530 mm	600 - 810 mm	370 - 520 mm	370 - 520 mm	740 mm