



## ASPIRE MAXI DAY CHAIR - UTILITY CHAIR

**\$1,401.00**



**SKU:** CHP208010

## PRODUCT DESCRIPTION

This sturdy bariatric seating support surface is an ideal bariatric seating solution with adjustable seat width, depth & backrest height.

## FEATURES

- Generous 310kg safe working load, a true bariatric solution.
- Ergonomic seat with anterior tilt to improve posture and reduce fatigue.
- Height, width and depth adjustability in seat and back for customised comfort.
- Easily manoeuvrable with push handle and transport wheels.
- Easy clean, breathable and comfortable vinyl upholstery and powder coated frame for effortless maintenance and longevity.

## SPECIFICATIONS

Name	Total Height	Safe Working Load	Material	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
Maxi Day Chair	1120 - 1440 mm	310 kg	Vinyl	45 kg	390 - 530 mm	600 - 810 mm	370 - 520 mm	370 - 520 mm	740 mm