



## **ASPIRE MAXI DAY CHAIR - UTILITY CHAIR**

\$1,277.00 Original price was:

\$1,277.00.<u>\$1,149.00</u>Current price is: \$1,149.00.



**SKU:** CHP208010

## **PRODUCT DESCRIPTION**

This sturdy bariatric seating support surface is an ideal bariatric seating solution with adjustable seat width, depth & backrest height.

## **FEATURES**

- Generous 310kg safe working load, a true bariatric solution.
- Ergonomic seat with anterior tilt to improve posture and reduce fatigue.
- Height, width and depth adjustability in seat and back for customised comfort.
- Easily manoeuvrable with push handle and transport wheels.
- Easy clean, breathable and comfortable vinyl upholstery and powder coated frame for effortless maintenance and longevity.

## **SPECIFICATIONS**

Name	Total Height	Safe Working Load	Material	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
	1120 - 1440	310 kg	Vinyl		390 - 530	600 - 810	370 - 520 mm	370 - 520	740 mm
Chair	mm				mm	mm		111111	