



ASPIRE MAXI SHOWER CHAIR



SKU: BTS118020

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aspire Maxi Shower Chair is a robust and comfortable solution designed for individuals who require additional support and stability during showering. This chair is particularly suited for older adults, those with mobility impairments, or anyone who may find standing in the shower challenging and unsafe.

FEATURES

- **Sturdy Construction:** Built with durable materials, the chair is designed to offer reliable support and withstand the wet conditions of a bathroom environment.
- **Comfortable Seating:** Features a wide, ergonomically designed seat that provides ample space and comfort, ensuring a safe and pleasant showering experience.
- **Adjustable Height:** The legs of the chair are adjustable, allowing the seat height to be modified to accommodate the needs and preferences of different users.
- **Non-Slip Rubber Feet:** Equipped with rubber feet to prevent slipping, ensuring the chair remains stable on wet bathroom floors.
- **Armrests and Backrest:** Comes with armrests and a backrest for additional support, aiding in maintaining balance and making it easier to rise from the chair.

The Aspire Maxi Shower Chair is an invaluable aid for those who face challenges with mobility or balance, offering a safe and secure way to enjoy showers independently. Its combination of practicality, durability, and comfort makes it a highly recommended choice for enhancing bathroom safety and independence.



SPECIFICATIONS

Code	Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load
BTS118020	Aspire MAXI Shower Chair	455 mm	575 - 775 mm	455 - 575 mm	310 kg

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773 brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777 elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au