



## ASPIRE NINJA CARBON WALKER

**\$889.00**



**SKU:** WAF705380-H

## PRODUCT DESCRIPTION

### Purpose & Benefits of the Aspire Ninja Carbon Walker

The Aspire Ninja Carbon Walker is a premium, lightweight mobility aid designed to offer superior manoeuvrability, comfort, and style. Made from durable carbon fibre, this walker is perfect for individuals seeking a sleek and practical solution for maintaining independence. Ideal for indoor and outdoor use, the Ninja Carbon Walker ensures stability and ease of use, providing reliable support for daily activities.

### Features

- **Ultra-Lightweight Carbon Fibre Frame:** The carbon fibre construction makes the walker exceptionally lightweight while maintaining durability and strength.
- **Height Adjustable Handles:** Easily adjustable handles ensure a tailored fit for each user, promoting proper posture and comfort.
- **Folding Design for Portability:** The compact folding mechanism allows for quick storage and transport, making it an excellent choice for travel or saving space at home.
- **Large Wheels for Smooth Navigation:** Equipped with large wheels to provide stability and a smooth ride on various surfaces, both indoors and outdoors.
- **Ergonomic Hand Grips with Brakes:** Features comfortable hand grips with integrated brakes for improved safety and control.
- **Comfortable Padded Seat and Backrest:** Offers a convenient place to rest during extended outings or daily errands.
- **Under-Seat Storage Bag:** Includes a secure, easily accessible storage bag to carry personal belongings or shopping items.

### Specifications

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Width	Seat To Floor Height	Seat To Footrest
WAF705380	Aspire Ninja	710 mm	625 mm	860 - 945 mm	150 kg	6.35 kg	460 mm	540 mm	540 mm
	Carbon Walker - Medium								
WAF705390	Aspire Ninja	710 mm	625 mm	860 - 1025 mm	150 kg	6.4 kg	460 mm	620 mm	620 mm
	Carbon Walker - Tall								

The Aspire Ninja Carbon Walker combines innovative materials with thoughtful design, providing a stylish and reliable mobility solution for users seeking lightweight support and convenience.