


ASPIRE NINJA CARBON WALKER



SKU: WAF705380-H

PRODUCT DESCRIPTION

Purpose & Benefits of the Aspire Ninja Carbon Walker

The Aspire Ninja Carbon Walker is a premium, lightweight mobility aid designed to offer superior manoeuvrability, comfort, and style. Made from durable carbon fibre, this walker is perfect for individuals seeking a sleek and practical solution for maintaining independence. Ideal for indoor and outdoor use, the Ninja Carbon Walker ensures stability and ease of use, providing reliable support for daily activities.

Features

- **Ultra-Lightweight Carbon Fibre Frame:** The carbon fibre construction makes the walker exceptionally lightweight while maintaining durability and strength.
- **Height Adjustable Handles:** Easily adjustable handles ensure a tailored fit for each user, promoting proper posture and comfort.
- **Folding Design for Portability:** The compact folding mechanism allows for quick storage and transport, making it an excellent choice for travel or saving space at home.
- **Large Wheels for Smooth Navigation:** Equipped with large wheels to provide stability and a smooth ride on various surfaces, both indoors and outdoors.
- **Ergonomic Hand Grips with Brakes:** Features comfortable hand grips with integrated brakes for improved safety and control.
- **Comfortable Padded Seat and Backrest:** Offers a convenient place to rest during extended outings or daily errands.
- **Under-Seat Storage Bag:** Includes a secure, easily accessible storage bag to carry personal belongings or shopping items.

Specifications

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Width	Seat To Floor Height	Seat To Footrest
WAF705380	Aspire Ninja Carbon Walker - Medium	710 mm	625 mm	860 - 945 mm	150 kg	6.35 kg	460 mm	540 mm	540 mm
WAF705390	Aspire Ninja Carbon Walker - Tall	710 mm	625 mm	860 - 1025 mm	150 kg	6.4 kg	460 mm	620 mm	620 mm

The Aspire Ninja Carbon Walker combines innovative materials with thoughtful design, providing a stylish and reliable mobility solution for users seeking lightweight support and convenience.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au