



## ASPIRE RAPHAEL QUATTRO POWERLIFT RECLINER CHAIR

**\$2,572.00 - \$3,195.00** Price range: \$2,572.00 through \$3,195.00



**SKU:** CHP2273-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Aspire Raphael Quattro Powerlift Recliner Chair is an embodiment of luxury and functionality, meticulously crafted to provide users with an unrivalled level of comfort and assistance. It stands out in the market for its robust construction and the unique Quattro motor system that allows for an infinite number of seating and reclining positions.

### FEATURES

- **Quattro Motor System:** Equipped with four independent motors, this chair offers users the flexibility to adjust the backrest, leg rest, lumbar, and head support independently, enabling them to find their perfect comfort position.
- **Powerlift Assistance:** The smooth and quiet lifting mechanism gently guides users from a sitting to a standing position, which is essential for those who have difficulty with transitions due to age or mobility restrictions.
- **Memory Settings:** With user-friendly controls, the chair can memorize the user's preferred positions, making it easy to return to the comfort setting of choice with just a touch of a button.
- **Ergonomic Design:** Every aspect of the chair is designed with ergonomics in mind, providing ample support to the back, neck, and legs, and reducing pressure points for long-term comfort.
- **Premium Upholstery:** The chair is typically covered with high-grade fabric or leather, offering both durability and a touch of elegance to any living space.
- **Integrated USB Port:** For added convenience, the chair often comes with a built-in USB port, allowing users to charge their devices without having to leave their comfortable spot.

The Aspire Raphael Quattro Powerlift Recliner Chair redefines what it means to relax. With its impeccable design choices, ranging from its advanced motor system to its plush upholstery options,

every detail is curated for the discerning user. If you're in pursuit of a chair that promises comfort, luxury, and longevity, look no further. This recliner chair is where innovation meets indulgence.

## SPECIFICATIONS

| Code      | Display Name                                   | Total Width | Total Height | Safe Working Load | Colour                 | Seat Depth | Seat Width | Seat To Floor Height | Seat To Footrest | Backrest Height |
|-----------|------------------------------------------------|-------------|--------------|-------------------|------------------------|------------|------------|----------------------|------------------|-----------------|
| CHP227290 | Aspire                                         |             |              |                   |                        |            |            |                      |                  |                 |
|           | Raphael                                        |             |              |                   |                        |            |            |                      |                  |                 |
|           | Quattro                                        |             |              |                   |                        |            |            |                      |                  |                 |
|           | Powerlift Recline Chair - Small - Bronze       | 740 mm      | 1060 mm      | 130 kg            | Bronze Microsuede      | 480 mm     | 480 mm     | 470 mm               | 470 mm           | 630 mm          |
| CHP227291 | Aspire                                         |             |              |                   |                        |            |            |                      |                  |                 |
|           | Raphael                                        |             |              |                   |                        |            |            |                      |                  |                 |
|           | Quattro                                        |             |              |                   |                        |            |            |                      |                  |                 |
|           | Powerlift Recline Chair - Small - Marina Grey  | 740 mm      | 1060 mm      | 130 kg            | Marina Grey Microsuede | 480 mm     | 480 mm     | 470 mm               | 470 mm           | 630 mm          |
| CHP227300 | Aspire                                         |             |              |                   |                        |            |            |                      |                  |                 |
|           | Raphael                                        |             |              |                   |                        |            |            |                      |                  |                 |
|           | Quattro                                        |             |              |                   |                        |            |            |                      |                  |                 |
|           | Powerlift Recline Chair - Medium - Bronze      | 780 mm      | 1090 mm      | 130 kg            | Bronze Microsuede      | 500 mm     | 520 mm     | 470 mm               | 470 mm           | 630 mm          |
| CHP227301 | Aspire                                         |             |              |                   |                        |            |            |                      |                  |                 |
|           | Raphael                                        |             |              |                   |                        |            |            |                      |                  |                 |
|           | Quattro                                        |             |              |                   |                        |            |            |                      |                  |                 |
|           | Powerlift Recline Chair - Medium - Marina Grey | 780 mm      | 1090 mm      | 130 kg            | Marina Grey Microsuede | 500 mm     | 520 mm     | 470 mm               | 470 mm           | 680 mm          |

| Code      | Display Name                                | Total Width | Total Height | Safe Working Load | Colour                 | Seat Depth | Seat Width | Seat To Floor Height | Seat To Footrest | Backrest Height |
|-----------|---------------------------------------------|-------------|--------------|-------------------|------------------------|------------|------------|----------------------|------------------|-----------------|
| CHP227302 | Aspire                                      |             |              |                   |                        |            |            |                      |                  |                 |
|           | Raphael                                     |             |              |                   |                        |            |            |                      |                  |                 |
|           | Quattro                                     |             |              |                   |                        |            |            |                      |                  |                 |
|           | Powerlift                                   |             |              |                   |                        |            |            |                      |                  |                 |
| CHP227319 | Recline Chair - Medium - Nordic Stone Vinyl | 780 mm      | 1090 mm      | 130 kg            | Nordic Stone Vinyl     | 500 mm     | 520 mm     | 470 mm               | 470 mm           | 680 mm          |
|           | Aspire                                      |             |              |                   |                        |            |            |                      |                  |                 |
|           | Raphael                                     |             |              |                   |                        |            |            |                      |                  |                 |
|           | Quattro                                     |             |              |                   |                        |            |            |                      |                  |                 |
| CHP227322 | Powerlift                                   | 930 mm      | 1050 mm      | 180 kg            | Bronze Microsuede      | 540 mm     | 650 mm     | 510 mm               | 510 mm           | 650 mm          |
|           | Recline Chair - Maxi - Bronze               |             |              |                   |                        |            |            |                      |                  |                 |
|           | Aspire                                      |             |              |                   |                        |            |            |                      |                  |                 |
|           | Raphael                                     |             |              |                   |                        |            |            |                      |                  |                 |
| CHP227322 | Quattro                                     |             |              |                   |                        |            |            |                      |                  |                 |
|           | Powerlift                                   | 930 mm      | 1050 mm      | 180 kg            | Marina Grey Microsuede | 540 mm     | 650 mm     | 510 mm               | 510 mm           | 650 mm          |
|           | Recline Chair - Maxi - Marina Grey          |             |              |                   |                        |            |            |                      |                  |                 |
|           | Aspire                                      |             |              |                   |                        |            |            |                      |                  |                 |