



ASPIRE RAPHAEL QUATTRO POWERLIFT RECLINER CHAIR

\$2,572.00 - \$3,195.00 Price range: \$2,572.00 through \$3,195.00



SKU: CHP2273-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aspire Raphael Quattro Powerlift Recliner Chair is an embodiment of luxury and functionality, meticulously crafted to provide users with an unrivalled level of comfort and assistance. It stands out in the market for its robust construction and the unique Quattro motor system that allows for an infinite number of seating and reclining positions.

FEATURES

- **Quattro Motor System:** Equipped with four independent motors, this chair offers users the flexibility to adjust the backrest, leg rest, lumbar, and head support independently, enabling them to find their perfect comfort position.
- **Powerlift Assistance:** The smooth and quiet lifting mechanism gently guides users from a sitting to a standing position, which is essential for those who have difficulty with transitions due to age or mobility restrictions.
- **Memory Settings:** With user-friendly controls, the chair can memorize the user's preferred positions, making it easy to return to the comfort setting of choice with just a touch of a button.
- **Ergonomic Design:** Every aspect of the chair is designed with ergonomics in mind, providing ample support to the back, neck, and legs, and reducing pressure points for long-term comfort.
- **Premium Upholstery:** The chair is typically covered with high-grade fabric or leather, offering both durability and a touch of elegance to any living space.
- **Integrated USB Port:** For added convenience, the chair often comes with a built-in USB port, allowing users to charge their devices without having to leave their comfortable spot.

The Aspire Raphael Quattro Powerlift Recliner Chair redefines what it means to relax. With its impeccable design choices, ranging from its advanced motor system to its plush upholstery options,

every detail is curated for the discerning user. If you're in pursuit of a chair that promises comfort, luxury, and longevity, look no further. This recliner chair is where innovation meets indulgence.

SPECIFICATIONS

Code	Display Name	Total Width	Total Height	Safe Working Load	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP227290	Aspire									
	Raphael									
	Quattro									
	Powerlift Recline Chair - Small - Bronze	740 mm	1060 mm	130 kg	Bronze Microsuede	480 mm	480 mm	470 mm	470 mm	630 mm
CHP227291	Aspire									
	Raphael									
	Quattro									
	Powerlift Recline Chair - Small - Marina Grey	740 mm	1060 mm	130 kg	Marina Grey Microsuede	480 mm	480 mm	470 mm	470 mm	630 mm
CHP227300	Aspire									
	Raphael									
	Quattro									
	Powerlift Recline Chair - Medium - Bronze	780 mm	1090 mm	130 kg	Bronze Microsuede	500 mm	520 mm	470 mm	470 mm	630 mm
CHP227301	Aspire									
	Raphael									
	Quattro									
	Powerlift Recline Chair - Medium - Marina Grey	780 mm	1090 mm	130 kg	Marina Grey Microsuede	500 mm	520 mm	470 mm	470 mm	680 mm

Code	Display Name	Total Width	Total Height	Safe Working Load	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP227302	Aspire									
	Raphael									
	Quattro									
	Powerlift									
CHP227319	Recline Chair - Medium - Nordic Stone Vinyl	780 mm	1090 mm	130 kg	Nordic Stone Vinyl	500 mm	520 mm	470 mm	470 mm	680 mm
	Aspire									
	Raphael									
	Quattro									
CHP227322	Powerlift Recline Chair - Maxi - Bronze	930 mm	1050 mm	180 kg	Bronze Microsuede	540 mm	650 mm	510 mm	510 mm	650 mm
	Aspire									
	Raphael									
	Quattro									
CHP227322	Powerlift Recline Chair - Maxi - Marina Grey	930 mm	1050 mm	180 kg	Marina Grey Microsuede	540 mm	650 mm	510 mm	510 mm	650 mm
	Aspire									
	Raphael									
	Quattro									