



ASPIRE REHAB RX ADVANCED TILT-IN-SPACE WHEELCHAIR

\$2,894.00 - \$3,315.00 Price range: \$2,894.00 through \$3,315.00



SKU: MWS4497-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Aspire Rehab RX Advanced Tilt-In-Space Wheelchair is a testament to modern mobility technology. Specifically designed for users with intricate requirements, this wheelchair guarantees optimum positioning, prolonged comfort, and superior skin management. Whether you need assistance with maintaining posture or are looking for an adaptive seating solution that addresses fatigue and skin health, the Rehab RX is a game-changer.

FEATURES

- **Positioning and Re-Positioning:** Built to cater to users with intricate needs, the Rehab RX ensures proper posture, optimises comfort, and significantly decreases user fatigue.
- **Tilt and Recline:** A unique feature that enables 26° of tilt and 30° of backrest recline. This not only helps in weight re-distribution but also preserves skin integrity, aids with head control, and augments overall user comfort.
- **Ergonomic Seating:** With a contoured, built-up backrest and seat cushion, the chair offers lateral support and promotes optimal midline positioning. The abductor and adductor build-ups at the seat cushion further enhance the user's seating posture.
- **Articulating Leg rests:** A design marvel that imitates natural knee movements during elevation. It ensures that the thigh remains in contact with the seating surface, providing unswerving support and maintaining skin contact with the cushion.
- **Breathable Design:** The multi-layered foam and breathable multi-stretch covers work hand in hand to redistribute pressure, boost comfort, and enhance skin microclimate.
- **Adaptable Design:** The wheelchair boasts adjustable seat depth and width, height-adjustable and removable armrests, and a versatile headrest, making it a truly customisable mobility

solution.

- **Safety Assured:** Crash Tested in accordance with AS/NZ3616.196 standards, this wheelchair is a fusion of reliability and functionality. The addition of anti-tippers as a standard feature further solidifies its stance on safety.

The Aspire Rehab RX Advanced Tilt-In-Space Wheelchair isn't just another mobility aid; it's a promise of enhanced quality of life. Tailored for users with complex requirements, it combines advanced technology with ergonomic design to deliver a wheelchair that's as comfortable as it is functional. From its pressure redistribution capabilities to its adaptive design, the Rehab RX is where mobility meets innovation.

SPECIFICATIONS

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Safe Working Load	Unit Weight	Colour	Material
MWS449630	Aspire Rehab RX Junior - 32cm - 37cm Wide	440 - 490 mm	32cm - 37cm (12.5" - 14.5")	460 - 710 mm	400 - 475 mm	180 kg	39.5 kg	Silver	Steel
MWS449710	Aspire Rehab RX - 41cm - 46cm Wide	440 - 530 mm	41cm - 46cm (16.1" - 18.1")	460 - 710 mm	400 - 475 mm	180 kg	41.0 kg	Silver	Steel
MWS449720	Aspire Rehab RX - 46cm - 51cm Wide	440 - 530 mm	46cm - 51cm (18.1" - 20")	460 - 710 mm	400 - 475 mm	180 kg	43.0 kg	Silver	Steel
MWS449740	Aspire Rehab RX - 51cm - 56cm Wide	470 - 560 mm	51cm - 56cm (20" - 22")	460 - 710 mm	400 - 475 mm	180 kg	47.2 kg	Silver	Steel