



ASPIRE REMBRANDT QUATTRO LIFT RECLINE CHAIR

\$6,054.00



SKU: CHP227550

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Welcome to the pinnacle of luxurious comfort – the Aspire Rembrandt Quattro Lift Recline Chair. This state-of-the-art chair isn't just a seating solution; it's an experience. Engineered for those who seek both luxury and utility, this chair offers a tailored comfort experience from head to toe. With features like a shiatsu massage, zero gravity positioning, and heating functions, it's more than just a recliner, it's a personal retreat.

FEATURES

- **Infinite Positioning:** Enjoy limitless seating options with independent control of the back, leg, lumbar, and headrest.
- **Shiatsu Roller Massage:** Indulge in a therapeutic shiatsu massage right from the comfort of your chair, providing stress relief and relaxation.
- **Memory Function:** Your comfort, remembered. Set and recall your favourite position for a consistent seating experience every time.
- **Integrated Lift and Recline:** Seamlessly transition from sitting to standing with the integrated lift and recline functions.
- **Extending Legrest:** Enhanced support for the lower legs ensures utmost comfort during relaxation.
- **Heating Function:** The cold will never bother you anyway with the chair's built-in heating system for the seat, lumbar, and backrest.
- **Built-in Conveniences:** Whether it's a cup holder for your favourite beverage or a wireless device charging spot integrated into the armrest, the Rembrandt Quattro has thought of everything.

- **Zero Gravity Positioning:** Experience the feeling of weightlessness while promoting pressure redistribution for optimum relaxation.
- **Power Backup:** Never be left without your chair's functionalities with the Lithium Ion Battery Backup, ensuring usage even during power outages or when moving between rooms.
- **Backlit Hand Control:** Modern controls featuring a backlit design with an integrated USB charging point, allowing you to keep your devices charged and within reach.

The Aspire Rembrandt Quattro Lift Recline Chair isn't just about sitting; it's about experiencing comfort at its peak. Every feature, from the massage to the heating functions, has been meticulously crafted to ensure users can relax, recharge, and rejuvenate. Whether you're watching a movie, reading a book, or just looking to escape into relaxation, the Rembrandt Quattro is your perfect companion.

SPECIFICATIONS

Code	Display Name	Total Width	Total Height	Safe Working Load	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height	User Height
CHP227550	Aspire Rembrandt Quattro Lift Recline Chair	1015 mm	1095 mm	159 Kg	Canyon Walnut	500 mm	560 mm	495 mm	495 mm	600 mm	1065 mm